



MEDMAIN

SURGICAL SUTURE

Global Medical Company, Search Medical



SEARCH MEDICAL

[61003] 16, Cheomdan venture so-ro 37beon-gil, Buk-gu, Gwangju, Korea

[970-73] Wolchul-dong, Buk-gu, Gwangju, Korea

Tel : +82-62-714-2699 Fax : +82-70-4855-2609 E-mail : sales@searchmedi.co.kr




MED MAIN

PRODUCTS Info

NYLON

Non-Absorbable Suture

Nylon Blue/Black, Mono filament



1.Raw material

Polyamide

2.Color

Blue, Black

3.Structure

Monofilament

4.Suture with Needle




[USP]

10/0, 9/0, 8/0, 7/0, 6/0, 5/0, 4/0, 3/0, 2/0, 1/0, 1, 2

[Length]

45cm, 75cm

5.Only Suture



[USP]

8/0, 7/0, 6/0, 5/0, 4/0, 3/0, 2/0, 1/0, 1, 2

[Length]

150cm, 50m, 75m, 100m

6.Validity


5years from date of manufacture

SILK

Non-Absorbable Suture

Silk White/Black, Braided

14p



1.Raw material

Silk

2.Color

Black, White

3.Structure

Multifilament

4.Suture with Needle



[USP]

8/0, 7/0, 6/0, 5/0, 4/0, 3/0, 2/0, 1/0, 1, 2

[Length]

45cm, 75cm

5.Only Suture



[USP]

8/0, 7/0, 6/0, 5/0, 4/0, 3/0, 2/0, 1/0, 1, 2, 3, 4

[Length]

150cm, 25m, 50m, 75m, 100m

6.Validity

5years from date of manufacture

SUPRACORE

Non-Absorbable Suture

Supracore White/Black, Multi filament

21p



1.Raw material

Supracore

5.Validity

5years from date of manufacture

2.Color

White, Black

3.Structure

Multi filament

4.Suture with Needle



[USP]

4/0

[Length]

45cm, 75cm

CODE INFORMATION

MEDMAIN CODE

REF Product code

NU20-HR10-45

1 2 3 4 5 6

01

1 봉합사 물질/ 색상

SW=Slik White

SB=Slik Black

NU=Nylon Blue

NB=Nylon Black


02

2 봉합사 두께 [USP]기준


03

3 봉합침 각도


F=5/8




H=1/2




Q=1/4



T=3/8



S=Straight



CODE INFORMATION

MEDMAIN CODE

REF Product code

NU20-HR10-45


1 2 3 4 5 6

04

4 봉합침 모양


R=Round body, Taper point

Round Body, Taper point ⊙




B= Round body, Blunt point

Round Body, Blunt Point ⊙




C=Reverse cutting

Reverse Cutting ▼




M=Micro point

Micro-point Reverse ▼



D=Diamond point

Diamond Point ⊗



S=Spatula

Spatula ▽



05

5 봉합침 길이 [mm]

06

6 봉합사 길이 [cm]

*봉합침 모양 다음“k”가 붙은 제품은 일반적인 봉합침보다 두꺼운 제품입니다.(NU20-HRK10-45)

MED MAIN



Polyamide

MEDMAIN With Surgical Suture Solution




Polyamide

LABELING ON THE PACKAGING



○ 상세설명

- | | |
|----------------------------|-------------------------|
| 1 제품명 | 2 봉합사 길이 |
| 3 제품성질 | 4 봉합사 사이즈(USP기준) |
| 5 제조사명 | 6 생산일 |
| 7 유효기간 | 8 제품코드 |
| 9 생산번호 | 10 봉합사 사이즈(EP기준) |
| 11 봉합침정보 (길이,각도,모양) | |







봉합침 코드 적용예시



- A=봉합침각도
- B=봉합침타입
- 1=봉합침길이













POLYAMIDE

Nylon Blue / Black, monofilament

Needle Product	Length [cm]	USP[Metric]				
		7/0 [0.5]	6/0 [0.7]	5/0 [1]	4/0 [1.5]	3/0 [2]
10mm Needle						
HR10 	45 Blue /Black	NU / NB 70 - HR10 - 45	NU / NB 60 - HR10 - 45			
2TR10 		NU / NB 70 - 2TR10 - 45	NU / NB 60 - 2TR10 - 45			
11mm Needle						
HC11 	45 Blue /Black		NU / NB 60 - HC11 - 45			
TCM11 		NU / NB 70 - TCM11 - 45	NU / NB 60 - TCM11 - 45			
2TCM11 		NU / NB 70 - 2TCM11 - 45	NU / NB 60 - 2TCM11 - 45			
12mm Needle						
TR12 	45 Blue /Black	NU / NB 70 - TR12 - 45	NU / NB 60 - TR12 - 45	NU / NB 50 - TR12 - 45		
TC12 		NU / NB 70 - TC12 - 45	NU / NB 60 - TC12 - 45	NU / NB 50 - TC12 - 45	NU / NB 40 - TC12 - 45	
2TC12 		NU / NB 70 - 2TC12 - 45	NU / NB 60 - 2TC12 - 45	NU / NB 50 - 2TC12 - 45	NU / NB 40 - 2TC12 - 45	
13mm Needle						
HR13 	45 Blue /Black		NU / NB 60 - HR13 - 45	NU / NB 50 - HR13 - 45		
TCM13 				NU / NB 50 - TCM13 - 45	NU / NB 40 - TCM13 - 45	
2TCM13 				NU / NB 50 - 2TCM13 - 45	NU / NB 40 - 2TCM13 - 45	











POLYAMIDE

Nylon Blue / Black, monofilament

Needle Product	Length [cm]	USP[Metric]				
		7/0 [0.5]	6/0 [0.7]	5/0 [1]	4/0 [1.5]	3/0 [2]
15mm Needle						
HR15 	45 Blue /Black			NU / NB 50 - HR15 - 45		
HC15 				NU / NB 50 - HC15 - 45	NU / NB 40 - HC15 - 45	NU / NB 30 - HC15 - 45
TR15 				NU / NB 50 - TR15 - 45		
16mm Needle						
TC16 	45 Blue /Black		NU / NB 60 - TC16 - 45	NU / NB 50 - TC16 - 45	NU / NB 40 - TC16 - 45	
17mm Needle						
HR17 	75 Blue /Black			NU / NB 50 - HR17 - 75	NU / NB 40 - HR17 - 75	NU / NB 30 - HR17 - 75
2HR17 				NU / NB 50 - 2HR17 - 75	NU / NB 40 - 2HR17 - 75	NU / NB 30 - 2HR17 - 75
18mm Needle						
2HR18 	45 Blue /Black	NU / NB 70 - 2HR18 - 45				
2HC18 		NU / NB 70 - 2HC18 - 45				
TR18 	75 Blue /Black			NU / NB 50 - TR18 - 75	NU / NB 40 - TR18 - 75	NU / NB 30 - TR18 - 75
2TR18 	45 Blue /Black	NU / NB 70 - 2TR18 - 45				
19mm Needle						
HR19 	45 Blue /Black				NU / NB 40 - HR19 - 45	NU / NB 30 - HR19 - 45
TC19 	45/75 Blue /Black			NU / NB 50 - TC19 - 45	NU / NB 40 - TC19 - 45/75	NU / NB 30 - TC19 - 45










POLYAMIDE

Nylon Blue / Black, monofilament

Needle Product	Length [cm]	USP[Metric]				
		7/0 [0.5]	6/0 [0.7]	5/0 [1]	4/0 [1.5]	3/0 [2]
20mm Needle						
 HR20	45/75 Blue /Black				NU / NB 40 - HR20 - 75	NU / NB 30 - HR20 - 45/75
 HC20	45 Blue /Black			NU / NB 50 - HC20 - 45	NU / NB 40 - HC20 - 45	NU / NB 30 - HC20 - 45
 2TR20		NU / NB 70 - 2TR20 - 45				
21mm Needle						
 TC21	45/75 Blue /Black				NU / NB 40 - TC21 - 45	NU / NB 30 - TC21 - 45/75
24mm Needle						
 HR24	45/75 Blue /Black	NU / NB 70 - HR24 - 45/75				NU / NB 30 - HR24 - 45/75
 TC24	45 Blue /Black	NU / NB 70 - TC24 - 45			NU/NB 40 - TC24 - 45	NU/NB 30 - TC24 - 45
Needle Product	Length [cm]	USP[Metric]				
		4/0 [1.5]	3/0 [2]	2/0 [3]	1/0 [4]	2 [5]
20mm Needle						
 HR20	75 Blue /Black		NU / NB 30 - HR20 - 45/75	NU / NB 20 - HR20 - 75		
23mm Needle						
 HC23	75 Blue /Black	NU / NB 40 - HC23 - 75	NU / NB 30 - HC23 - 75	NU / NB 20 - HC23 - 75		
24mm Needle						
 HR24	45/75 Blue /Black		NU / NB 30 - HR24 - 45/75	NU / NB 20 - HR24 - 75		
 TC24	45 Blue /Black	NU / NB 40 - TC24 - 45	NU / NB 30 - TC24 - 45	NU / NB 20 - TC24 - 45		





POLYAMIDE

Nylon Blue / Black, monofilament

Needle Product	Length [cm]	USP[Metric]				
		4/0 [1.5]	3/0 [2]	2/0 [3]	1/0 [3.5]	1 [4]
26mm Needle						
 HR26	75 Blue /Black	NU / NB 40 - HR26 - 75	NU / NB 30 - HR26 - 75	NU / NB 20 - HR26 - 75	NU / NB 10 - HR26 - 75	
 2HR26		NU / NB 40 - 2HR26 - 75	NU / NB 30 - 2HR26 - 75	NU / NB 20 - 2HR26 - 75		
 HC26			NU / NB 30 - HC26 - 75	NU / NB 20 - HC26 - 75		
 TR26		NU / NB 40 - TR26 - 75	NU / NB 30 - TR26 - 75	NU / NB 20 - TR26 - 75		
30mm Needle						
 HR30	75 Blue /Black		NU / NB 30 - HR30 - 75	NU / NB 20 - HR30 - 75	NU / NB 10 - HR30 - 75	NU / NB 01 - HR30 - 75
 TC30			NU / NB 30 - TC30 - 75	NU / NB 20 - TC30 - 75	NU / NB 10 - TC30 - 75	NU / NB 01 - TC30 - 75
32mm Needle						
 TR32	75 Blue /Black			NU / NB 20 - TR32 - 75	NU / NB 10 - TR32 - 75	NU / NB 01 - TR32 - 75
37mm Needle						
 HR37	75 Blue /Black		NU / NB 30 - HR37 - 75	NU / NB 20 - HR37 - 75	NU / NB 10 - HR37 - 75	NU / NB 01 - HR37 - 75
39mm Needle						
 TC39	45/75 Blue /Black			NU / NB 20 - TC39 - 45/75	NU / NB 10 - TC39 - 75	NU / NB 01 - TC39 - 75

POLYAMIDE

Nylon Blue / Black, monofilament

Needle Product	Length [cm]	USP[Metric]				
		3/0 [2]	2/0 [3]	1/0 [3.5]	1 [4]	2 [5]
40mm Needle						
 HR40	75/100 Blue /Black		NU / NB 20 - HR40-75/100			
45mm Needle						
 TR45	100 Blue /Black		NU / NB 20 - TR45-100			
48mm Needle						
 HR48	75 Blue /Black		NU / NB 20 - HR48 - 75	NU / NB 10 - HR48 - 75	NU / NB 01 - HR48 - 75	NU / NB 02 - HR48 - 75
60mm Needle						
 SC60	75 Blue /Black		NU / NB 20 - SC60 - 75	NU / NB 10 - SC60 - 75	NU / NB 01 - SC60 - 75	

Without Needle	Color	USP[Metric]										
		8/0 [0.4]	7/0 [0.5]	6/0 [0.7]	5/0 [1]	4/0 [1]	3/0 [2]	2/0 [3]	1/0 [3.5]	1 [3.5]	2 [5]	
Without needle suture												
150cm	Blue /Black	NU / NB 80 - 150CM	NU / NB 70 - 150CM	NU / NB 60 - 150CM	NU / NB 50 - 150CM	NU / NB 40 - 150CM	NU / NB 30 - 150CM	NU / NB 20 - 150CM	NU / NB 10 - 150CM	NU / NB 01 - 150CM	NU / NB 02 - 150CM	
50m											NU / NB 02 - 50M	
75m										NU / NB 10 - 75M	NU / NB 01 - 75M	NU / NB 02 - 75M
100m			NU / NB 80 - 100M	NU / NB 70 - 100M	NU / NB 60 - 100M	NU / NB 50 - 100M	NU / NB 40 - 100M	NU / NB 30 - 100M	NU / NB 20 - 100M	NU / NB 10 - 100M	NU / NB 01 - 100M	

POLYAMIDE

MEDMAIN


MEDMAIN With Surgical Suture Solution



MEDMAIN With Surgical Suture Solution


SILK

LABELING ON THE PACKAGING



○ 상세설명

- | | |
|----------------------------|-------------------------|
| 1 제품명 | 2 봉합사 길이 |
| 3 제품성질 | 4 봉합사 사이즈(USP기준) |
| 5 제조사명 | 6 생산일 |
| 7 유효기간 | 8 제품코드 |
| 9 생산번호 | 10 봉합사 사이즈(EP기준) |
| 11 봉합침정보 (길이,각도,모양) | |
















봉합침 코드 적용예시

A B 1 1











A=봉합침각도

B=봉합침타입

1=봉합침길이











Needle Product	Length [cm]	USP[Metric]					
		8/0 [0.4]	7/0 [0.5]	6/0 [0.7]	5/0 [1]	4/0 [1.5]	
10mm Needle							
HR10 	45 White /Black		SW / SB 70 - HR10 - 45	SW / SB 60 - HR10 - 45			
TR10 			SW / SB 70 - TR10 - 45	SW / SB 60 - TR10 - 45			
2TR10 			SW / SB 70 - 2TR10 - 45	SW / SB 60 - 2TR10 - 45			
TC10 	45/75 White /Black				SW / SB 50 - TC10 - 45		
11mm Needle							
TCM11 	45 White /Black		SW / SB 70 - TCM11 - 45	SW / SB 60 - TCM11 - 45			
2TCM11 				SW / SB 70 - 2TCM11 - 45	SW / SB 60 - 2TCM11 - 45		
12mm Needle							
TR12 	45 White /Black		SW / SB 70 - TR12 - 45	SW / SB 60 - TR12 - 45	SW / SB 50 - TR12 - 45		
2TR12 				SW / SB 70 - 2TR12 - 45	SW / SB 60 - 2TR12 - 45	SW / SB 50 - 2TR12 - 45	
TC12 				SW / SB 70 - TC12 - 45	SW / SB 60 - TC12 - 45	SW / SB 50 - TC12 - 45	SW / SB 40 - TC12 - 45
2TC12 				SW / SB 70 - 2TC12 - 45	SW / SB 60 - 2TC12 - 45	SW / SB 50 - 2TC12 - 45	SW / SB 40 - 2TC12 - 45
13mm Needle							
HR13 	45 White /Black			SW / SB 60 - HR13 - 45	SW / SB 50 - HR13 - 45		
TCM13 					SW / SB 50 - TCM13 - 45	SW / SB 40 - TCM13 - 45	
2TCM13 					SW / SB 50 - 2TCM13 - 45	SW / SB 40 - 2TCM13 - 45	

Silk whith/black, Braided






Needle Product	Length [cm]	USP[Metric]				
		7/0 [0.5]	6/0 [0.7]	5/0 [1]	4/0 [1.5]	3/0 [2]
15mm Needle						
 HR15	75 White /Black				SW / SB 40 - HR15 - 75	
 HC15	45 White /Black			SW / SB 50 - HC15 - 45	SW / SB 40 - HC15 - 45	SW / SB 30 - HC15 - 45
 TR15	45/75 White /Black				SW / SB 40 - TR15 - 45 / 75	
16mm Needle						
 TC16	45 White /Black		SW / SB 60 - TC16 - 45	SW / SB 50 - TC16 - 45	SW / SB 40 - TC16 - 45	
17mm Needle						
 HR17	75 White /Black			SW / SB 50 - HR17 - 75	SW / SB 40 - HR17 - 75	SW / SB 30 - HR17 - 75
 2HR17				SW / SB 50 - 2HR17 - 75	SW / SB 40 - 2HR17 - 75	SW / SB 30 - 2HR17 - 75
18mm Needle						
 TR18	75 White /Black			SW / SB 50 - TR18 - 75	SW / SB 40 - TR18 - 75	SW / SB 30 - TR18 - 75
19mm Needle						
 TC19	45 White /Black			SW / SB 50 - TC19 - 45	SW / SB 40 - TC19 - 45	SW / SB 30 - TC19 - 45
Needle Product	Length [cm]	USP[Metric]				
		5/0 [1]	4/0 [1.5]	3/0 [2]	2/0 [3]	1/0 [3.5]
20mm Needle						
 HR20	45/75 White /Black		SW / SB 40 - HR20 - 45 / 75	SW / SB 30 - HR20 - 45 / 75	SW / SB 20 - HR20 - 45 / 75	
 HC20	45 White /Black	SW / SB 50 - HC20 - 45	SW / SB 40 - HC20 - 45	SW / SB 30 - HC20 - 45		

SILK

Silk whith/black, Braided

Needle Product	Length [cm]	USP[Metric]				
		5/0 [1]	4/0 [1.5]	3/0 [2]	2/0 [3]	1/0 [3.5]
21mm Needle						
 TR21	75 White /Black			SW / SB 30 - TR21 - 75		
 TC21	45 White /Black	SW / SB 50 - TC21 - 45	SW / SB 40 - TC21 - 45	SW / SB 20 - TC21 - 45		
23mm Needle						
 HC23	75 White /Black		SW / SB 40 - HC23 - 75	SW / SB 30 - HC23 - 75	SW / SB 20 - HC23 - 75	
24mm Needle						
 TC24	45 White /Black		SW / SB 40 - TC24 - 45	SW / SB 30 - TC24 - 45	SW / SB 20 - TC24 - 45	
26mm Needle						
 HR26	45/75 White /Black		SW / SB 40 - HR26 - 45 / 75	SW / SB 30 - HR26 - 45 / 75	SW / SB 20 - HR26 - 45 / 75	SW / SB 10 - HR26 - 45 / 75
 2HR26	75 White /Black		SW / SB 40 - 2HR26 - 75	SW / SB 30 - 2HR26 - 75	SW / SB 20 - 2HR26 - 75	
 HC26				SW / SB 30 - HC26 - 75	SW / SB 20 - HC26 - 75	
 TR26			SW / SB 40 - TR26 - 75	SW / SB 30 - TR26 - 75	SW / SB 20 - TR26 - 75	
30mm Needle						
 HR30	75 White /Black	SW / SB 30 - HR30 - 75	SW / SB 20 - HR30 - 75	SW / SB 10 - HR30 - 75	SW / SB 01 - HR30 - 75	
 TC30	45/75 White /Black	SW / SB 30 - TC30 - 75	SW / SB 20 - TC30 - 75	SW / SB 10 - TC30 - 75	SW / SB 01 - TC30 - 45	


Silk whith/black, Braided

Needle Product	Length [cm]	USP[Metric]												
		3/0 [2]	2/0 [3]	1/0 [3.5]	1 [4]	2 [5]								
32mm Needle														
 TR32	75 White /Black		SW / SB 20 - TR32 - 75	SW / SB 10 - TR32 - 75	SW / SB 01 - TR32 - 75									
37mm Needle														
 HR37	75 White /Black	SW / SB 30 - HR37 - 75	SW / SB 20 - HR37 - 75	SW / SB 10 - HR37 - 75	SW / SB 01 - HR37 - 75									
39mm Needle														
 TC39	75 White /Black		SW / SB 20 - TC39 - 75	SW / SB 10 - TC39 - 75	SW / SB 01 - TC39 - 75									
48mm Needle														
 HR48	75 White /Black		SW / SB 20 - HR48 - 75	SW / SB 10 - HR48 - 75	SW / SB 01 - HR48 - 75	SW / SB 02 - HR48 - 75								
60mm Needle														
 SC60	75 White /Black		SW / SB 20 - SC60 - 75	SW / SB 10 - SC60 - 75	SW / SB 10 - SC60 - 75									
Without Needle	Color	USP[Metric]												
		8/0 [0.4]	7/0 [0.5]	6/0 [0.7]	5/0 [1]	4/0 [1]	3/0 [2]	2/0 [3]	1/0 [3.5]	1 [3.5]	2 [5]	3 [6]	4 [6]	
Without needle suture														
150cm	White /Black	SW / SB 80 - 150CM	SW / SB 70 - 150CM	SW / SB 60 - 150CM	SW / SB 50 - 150CM	SW / SB 40 - 150CM	SW / SB 30 - 150CM	SW / SB 20 - 150CM	SW / SB 10 - 150CM	SW / SB 01 - 150CM	SW / SB 02 - 150CM	SW / SB 03 - 150CM	SW / SB 04 - 150CM	
													SW / SB 04 - 25M	
												SW / SB 02 - 50M	SW / SB 03 - 50M	
										SW / SB 10 - 75M	SW / SB 01 - 75M			
		SW / SB 80 - 100M	SW / SB 70 - 100M	SW / SB 60 - 100M	SW / SB 50 - 100M	SW / SB 40 - 100M	SW / SB 30 - 100M	SW / SB 20 - 100M	SW / SB 10 - 100M					


SILK
MEDMAIN
MEDMAIN With Surgical Suture Solution

Supracore

MEDMAIN With Surgical Suture Solution



SUPRACORE LABELING ON THE PACKAGING



○ 상세설명

1 제품명

3 제품성질

5 제조사명

7 유효기간

9 생산번호

11 봉합침정보 (길이,각도,모양)

2 봉합사 길이

4 봉합사 사이즈(USP기준)

6 생산일

8 제품코드

10 봉합사 사이즈(EP기준)

3/8 circle
Reverse Cutting

→ 봉합침 설명

TC19

→ 봉합침 코드



→ 봉합침 타입

19mm

→ 봉합침 길이



→ 봉합침 모양

봉합침 코드 적용예시

A B 1 1



A=봉합침각도

B=봉합침타입

1=봉합침길이

SUPRACORE

Supracore whith/black, Multi filament

Needle Product	Length [cm]	USP[Metric]
		4/0 [1.5]
19mm Needle		
 TC19	45 Black	AB40-TC19-45
	45 White	AB40-TC19-75
	75 Black	AW40-TC19-45
	75 White	AW40-TC19-75
 TCK19	45 Black	AB40-TCK19-45
	45 White	AW40-TCK19-45
	75 Black	AB40-TCK19-75
	75 White	AW40-TCK19-75

SUPRACORE

MEDMAIN

MEDMAIN With Surgical Suture Solution

MEDMAIN

With Surgical Suture Solution



SEARCH MEDICAL

[61003] 16, Cheomdan venture so-ro 37beon-gil, Buk-gu, Gwangju, Korea
[970-73] Wolchul-dong, Buk-gu, Gwangju, Korea